



SVANOVIR® BRSV-Ab

Detecting antibodies specific to Bovine Respiratory Syncytial Virus (BRSV) in cattle

SUMMARY | SVANOVIR® BRSV-Ab detects antibodies specific to BRSV in serum/plasma and milk from individuals, and in bulk milk samples. The assay provides accurate results, which are essential for the evaluation of maternal immunity in calves before vaccination, for the identification of cows exposed to BRSV, and herds free from infection.



YOUR CHALLENGE is a virus that paves the way for bacterial infection

BRSV is ubiquitous and spreads rapidly in dense cattle populations, especially those with low herd immunity. Secondary opportunistic bacterial or mycoplasmal infections follow exposure to BRSV and frequently aggravate the clinical symptoms. BRSV disease can be often severe in calves after a decline of their maternal antibodies and in previously non-exposed animals.

YOUR GOAL is to detect the virus behind the respiratory disease

There are several different viruses that can cause respiratory disease. Besides BRSV, those of major significance include Bovine Herpesvirus 1, Parainfluenza 3 and Bovine Coronavirus. Farmers face economic loss owing to decreased performance and considerable consumption of medication to prevent secondary bacterial infection. The most important factor in dealing with respiratory diseases is changing management routines to improve herd health.

High specificity – reliable results for controlling the spread of BRSV

Effective screening of bulk tank milk

Enables informed decision making in management routines

Screening of calves – determine the optimal time for vaccination after decline of maternal antibodies

ASSAY OVERVIEW



SVANOVIR® BRSV-Ab

Species	Bovine		
Samples	Serum/plasma Milk, individual and bulk milk		
Type	Indirect ELISA detecting IgG ₁		
Article number	Samples*	Plates	Format
104888	88	2	Strips

* Samples: Max. number of samples for analysis, wells for kit controls excluded.

SVANOVIR® BRSV-Ab is a valuable tool for evaluating exposure to BRSV in cattle and provides essential information for management routines and movement of livestock.

Dual functionality – the same kit can be used for analysis of serum/plasma and milk samples

Work effective – single dilution assay


Easy protocol and ready-to-use reagents

High quality – validated and manufactured under strict ISO 9001:2008 standardised procedures in Sweden

Multilingual labels

YOUR SUPPORT

From 9am-16pm CET call:

 +46 18 65 49 15

 customer.service@svanova.com

Boehringer Ingelheim Svanova
Box 1545
SE-751 45 Uppsala, Sweden

www.svanova.com

svanova 

PERFORMANCE CHARACTERISTICS SVANOVIR® BRSV-Ab

On a large number of serum samples from two different populations (herds free of infection and herds with BRSV infection) in Sweden and England, the assay detected antibodies to BRSV accurately. Good agreement is seen between serum and milk antibody levels. Additionally the assay detects differences in antibody levels (3-fold increase of OD values) in paired diagnostic samples (Elvander et al., 1995). SVANOVIR® BRSV-Ab has also been successfully used for monitoring infection of BRSV in calves (Hägglund et al., 2006).

Serum from dairy cows
n = 157^a

Test	Reference ELISA		
	Positives	Negatives	Total
SVANOVIR® BRSV-Ab			
Positives	93	0	93
Negatives	6	58	64
Total	99	58	157
	Sensitivity	94%	
	Specificity	100%	

Samples from (a) Sweden and England from populations with history of respiratory disease

References

Elvander M., Edwards S, Näslund IS, Linde N. (1995): Evaluation and Application of an Indirect ELISA for the Detection of Antibodies to Bovine Respiratory Syncytial Virus in Milk, Bulk Milk, and Serum. J VET Diagn Invest 1995 7: 177

COMPLIMENTARY DIAGNOSTIC SOLUTIONS

SVANOVIR®BCV-Ab - Detecting antibodies specific to Bovine Coronavirus

SVANOVIR®PIV3-Ab - Detecting antibodies specific to Bovine Parainfluenza Virus Type 3

SVANOVIR®IBR-Ab - Detecting antibodies specific to Bovine Herpesvirus 1