



## SVANOVIR® PRV gE-Ab

# Controlling Pseudorabies / Aujeszky's disease in vaccinated swine populations

**SUMMARY** | SVANOVIR® PRV gE-Ab is a blocking ELISA for the detection of Pseudorabies / Aujeszky's infection in vaccinated swine populations. The high specificity enables discrimination of serological response to gE-deleted vaccine strains from that of field virus. This is a great benefit for use in eradication as well as control programmes and for certifying swine for import and export.



### **YOUR CHALLENGE** is a persisting herpes virus

The severity of the clinical manifestation in swine is age dependant, where severe disease with fatal outcome is more often seen in young piglets. In adults the disease is fairly mild and after clinical recovery persistently infected animals are a major risk for virus transmission. Control of this globally occurring disease is achieved by stamping out infected animals and/or by vaccinating swine populations at risk.

### **YOUR GOAL** is to discriminate between antibodies from vaccination vs. field infection in animals

Efficient vaccines reduce or prevent clinical signs without necessarily preventing virus replication in individuals and herds. Marker vaccines are part of the DIVA concept (differentiating infected from vaccinated individuals) and specific diagnostic tests are an essential complement to discriminate a serological response induced by vaccination from those of field infection.



### **ASSAY OVERVIEW**

#### SVANOVIR® PRV gE-Ab

<b>Species</b>	Porcine		
<b>Samples</b>	Serum/plasma, whole blood (filter paper)*		
<b>Type</b>	Blocking ELISA based on full antigen, anti-gE(gI) monoclonal antibodies		
<b>Article number</b>	<b>Samples**</b>	<b>Plates</b>	<b>Format</b>
104902	184	2	Strips
104904	920	10	Plates

\* Extra protocol available on request

\*\* Samples: Max number of samples for analysis, wells for kit controls excluded

### **Reliable test results**


**High performing test** with high sensitivity and specificity

**Discriminates vaccinated animals** from naturally infected

**Detects carrier animals** from naturally infected

**Standardised** against the reference serum ADV-1

**Field proven** and used in eradication and control programmes in Europe



The SVANOVIR® PRV gE-Ab assay has demonstrated excellent performance and is a field-proven tool in the control of Aujeszky's disease in vaccinated swine populations.


**Effective handling and convenient in use** with flexible protocol, ready-to-use reagents and alternative incubation time

**High quality** thoroughly validated and manufactured under strict ISO 9001:2008 standardised procedures in Sweden

**Multilingual kit insert and labels**

## YOUR SUPPORT

From 9am-4pm CET call:

 **+46 18 65 49 15**

 **customer.service@svanova.com**

Boehringer Ingelheim Svanova  
Box 1545  
SE-751 45 Uppsala, Sweden

**www.svanova.com**

svanova 

## PERFORMANCE CHARACTERISTICS SVANOVIR® PRV gE-Ab

In several studies including well-defined samples from different epidemiological subpopulations SVANOVIR® PRV gE-Ab assay has demonstrated superior performance. The test discriminates between antibodies from vaccination and field infection. False positive results is minimised due to the high specificity of the test. The high sensitivity allows correct detection of positive animals.

Samples	Sensitivity	Specificity	Reference method
Commercial pigs n= 1392 <sup>a</sup>	100%	99.6%	VNT
Naturally infected, vaccinated, Non-vaccinated herds n= 999 <sup>b</sup>	98.9%*	99.6%	DANISH BLOCKING ELISA**
Well defined samples from National re- ference Center <sup>c</sup>	97%	100%	In-house ELISA IZSLER Brescia
Commercial pigs n= 155 <sup>d</sup>	n.a%	100%	VNT

Samples originating from <sup>a</sup>Sweden, Yugoslavia, Germany, <sup>b</sup>Denmark, <sup>c</sup>Italy- National Reference Center for Aujeszky Disease - (IZLER), Brescia, <sup>d</sup>Germany\*, Virus neutralization test, \*\*Sørensen and Lei (1986), J. Virol. Methods

In a study performed by the National Reference Centre for Aujeszky's disease in Italy - on 46 well-defined samples - SVANOVIR® PRV gE-Ab showed excellent agreement with their reference method (in-house ELISA).

**References:** Sørensen KJ, Lei JC: 1986, Aujeszky's disease: blocking ELISA for the detection of serum antibodies. J Virol Methods 13:171-181

## Complementary product

<b>SVANOVIR® PRV gB-Ab</b>	Controlling Pseudorabies/Aujeszky's disease in non-vaccinated swine populations
----------------------------	---