


# **PRODUCT CATALOGUE**

## Content

Bovine	5
Ovine	9
Porcine	11
Equine	13
Avian	14







## SBV

### Schmallenberg Virus Antibody Test

An indirect ELISA for the detection of specific antibodies to Schmallenberg Virus in bovine serum, plasma and milk samples (individual and bulk tank milk).

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® SBV-Ab	113124	strips	2	96	88

## IBR

### Infectious Bovine Rhinotracheitis Antibody Test

An indirect ELISA for the detection of specific antibodies to Infectious Bovine Rhinotracheitis in bovine serum, plasma and milk samples (individual and bulk tank milk).

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® IBR-Ab	104878	strips	2	96	88
SVANOVIR® IBR-Ab	104879	plate	10	480	440
SVANOVIR® IBR-Ab (screen)	104880	strips	10	960	920

## BVDV

### Bovine Virus Diarrhoea Virus Antibody Test

An indirect ELISA for the detection of specific antibodies to Bovine Virus Diarrhoea Virus in bovine serum, plasma and milk samples (individual and bulk tank milk).

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® BVDV-Ab	104881	strips	2	96	88
SVANOVIR® BVDV-Ab	104882	plate	10	480	440
SVANOVIR® BVDV-Ab (screen)	104883	strips	10	960	920

### Bovine Virus Diarrhoea Virus p80 Antibody Test

A blocking ELISA for the detection of specific antibodies to Bovine Virus Diarrhoea p80 (NS3) in bovine, caprine and ovine serum, plasma and milk samples (individual and bulk tank milk).

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® BVDV p80-Ab	104911	strips	2	192	184

## BLV

### Bovine Leukaemia Virus gp51 Antibody Test

An indirect ELISA for the detection of specific antibodies to Bovine Leukaemia Virus gp51 in bovine serum, plasma and milk samples (individual and bulk tank milk).

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® BLV gp51-Ab	104884	strips	2	96	88
SVANOVIR® BLV gp51-Ab	104885	plate	10	480	440
SVANOVIR® BLV gp51-Ab (screen)	104886	strips	10	960	920

## BCV

### Bovine Coronavirus Antibody Test

An indirect ELISA for the detection of specific antibodies to Bovine Coronavirus in bovine serum and milk samples (individual and bulk tank milk).

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® BCV-Ab	104887	strips	2	96	88

## BRSV

### **Bovine Respiratory Syncytial Virus Antibody Test**

An indirect ELISA for the detection of specific antibodies to Bovine Respiratory Syncytial Virus in bovine serum and milk samples (individual and bulk tank milk).

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® BRSV-Ab	104888	strips	2	96	88

## PIV3

### **Bovine Parainfluenza Virus Type 3 Antibody Test**

An indirect ELISA for the detection of specific antibodies to Bovine Parainfluenza Virus Type 3 in bovine serum and milk samples.

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® PIV3-Ab	104889	strips	2	96	88

## BRUCELLA

### **Brucella I-ELISA Antibody Test**

An indirect ELISA for the detection of specific antibodies to *Brucella* in bovine serum, plasma and milk samples (individual and bulk tank milk).

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® <i>Brucella</i> -Ab I-ELISA	104890	strips	10	960	920

### **Brucella C-ELISA Antibody Test**

A competitive ELISA for the detection of specific antibodies to *Brucella* in bovine, ovine, caprine and porcine serum and plasma samples; which is capable of discriminating between naturally infected and S19 vaccinated cattle.

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® <i>Brucella</i> -Ab C-ELISA	104893	strips	2	192	176
SVANOVIR® <i>Brucella</i> -Ab C-ELISA	104894	plate	10	960	880

## FMDV

### **Foot and Mouth Disease Virus 3ABC Antibody Test – Ruminant**

An indirect ELISA for the detection of specific antibodies to Foot and Mouth Disease Virus 3ABC in bovine, ovine and caprine serum samples; which is capable of discrimination between naturally infected and vaccinated animals.

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® FMDV 3ABC-Ab ruminant	104895	plate	10	480	440

### **IAEA Standard Bovine FMDV Positive Control Sera**

A panel of reference sera for FMDV assays (lyophilised).

Product name	Article No.	Quantity
IAEA standard FMDV type A	113126	0.5 ml
IAEA standard FMVDV type O	113128	0.5 ml
IAEA standard FMDV type Asia1	113129	0.5 ml
IAEA standard FMDV type SAT1	113130	0.5 ml
IAEA standard FMDV type SAT2	113131	0.5 ml
IAEA standard FMDV type SAT3	113132	0.5 ml

Tests - number of tests. Samples - number of samples, wells for kit controls excluded.

## FMDV

### Foot and Mouth Disease Virus Antigen Test

A penside test for the detection of Foot and Mouth Disease Virus in swab and tissue samples.

Product name	Article No.	No of tests
SVANODIP® FMDV-Ag	104737	20 test devices

### Foot and Mouth Disease Virus Tissue Sample Extraction Kit

A tissue sample extraction kit designed to extract the antigen from tissues.

Product name	Article No.	No of tests
SVANODIP® FMDV-Ag Extraction	104872	20 extraction vials

## RPV

### Rinderpest Virus Antigen Test

A penside test for the detection of Rinderpest virus in eye swab samples.

Product name	Article No.	No of tests
SVANODIP® RPV-Ag	104738	50 test devices

## FASCIOLA

### Fasciola hepatica Antibody Test

An indirect ELISA for the detection of antibodies to *Fasciola hepatica* in bovine serum, plasma, meat juice and milk samples (individual and bulk tank milk), and for the classification of herds and individual animals predicting the economic consequences of infection.

Product name	Article no.	Format	Plates	Tests	Samples
SVANOVIR® <i>F. hepatica</i> -Ab	104896	strips	2	192	184

## OSTERTAGIA

### Ostertagia ostertagi Antibody Test

An indirect ELISA for the detection of specific antibodies to *Ostertagia ostertagi* in bovine bulk tank milk.

Product name	Article no.	Format	Plates	Tests	Samples
SVANOVIR® <i>O. ostertagi</i> -Ab	104897	strips	2	192	184

## NEOSPORAS

### Neospora caninum iscom ELISA Antibody Test

An indirect ELISA for the detection of specific antibodies to *Neospora caninum* in bovine serum and milk samples (individual and bulk tank milk).

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® <i>Neospora</i> -Ab	104898	strips	2	192	184

## ANAPLASMA

### Anaplasma marginale Antibody Test

An indirect ELISA for the detection of specific antibodies to *Anaplasma marginale* in bovine serum samples.

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® <i>A. marginale</i> -Ab	104899	strips	2	192	184

## BABESIA

### Babesia bigemina Antibody Test

An indirect ELISA for the detection of specific antibodies to *Babesia bigemina* in bovine serum samples.

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® <i>B. bigemina</i> -Ab	104900	strips	2	192	184





# OVINE / CAPRINE

## SBV

### Schmallenberg Virus Antibody Test

An indirect ELISA for the detection of specific antibodies to Schmallenberg Virus in ovine and caprine serum and plasma samples.

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® SBV-Ab	113124	strips	2	96	88

## BDV

### Border Disease Virus Antibody Test

An indirect ELISA for the detection of antibodies to Border Disease Virus in ovine and caprine serum and milk samples.

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® BDV-Ab	104910	strips	2	96	88

## BRUCELLA

### Brucella I-ELISA Antibody Test (small ruminants)

An indirect ELISA for the detection of specific antibodies to *Brucella* in ovine and caprine serum and plasma samples.

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® <i>Brucella</i> -Ab I-ELISA o/c	104892	strips	2	192	184

### Brucella C-ELISA Antibody Test

An competitive ELISA for the detection of specific antibodies to *Brucella* in ovine, caprine, bovine and porcine serum and plasma samples.

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® <i>Brucella</i> -Ab C-ELISA	104893	strips	2	192	176
SVANOVIR® <i>Brucella</i> -Ab C-ELISA	104894	plate	10	960	880

## FMDV

### Foot and Mouth Disease Virus 3ABC Antibody Test – Ruminant

An indirect ELISA for the detection of specific antibodies to Foot and Mouth Disease Virus 3ABC in bovine, ovine and caprine serum samples; which is capable of discrimination between naturally infected and vaccinated animals.

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® FMDV 3ABC-Ab ruminant	104895	plate	10	480	440



## LAWSONIA

### Porcine *Lawsonia intracellularis* Antibody Test

A blocking ELISA for the detection of specific antibodies to *Lawsonia intracellularis* in porcine serum and plasma samples.

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® <i>L. intracellularis</i> /Ileitis-Ab*	122275	strips	2	192	184

\*SVANOVIR® *L. intracellularis*/Ileitis-Ab is the compliant successor of the bioScreen Ileitis Antibody ELISA.

## TGEV/PRCV

### Transmissible Gastroenteritis Virus/Porcine Respiratory Coronavirus Differentiating test

A blocking ELISA for the detection and/or differentiation of specific antibodies to Transmissible Gastroenteritis and Porcine Respiratory Coronavirus in porcine serum samples.

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® TGEV/PRCV-Ab	104909	strips	2	96	88

## SALMONELLA

### Porcine *Salmonella* Antibody Test

An indirect ELISA for the detection of specific antibodies to *Salmonella* in pig serum and meat juice samples.

Product name	Article No.	Format	Plates	Tests	Samples
VETSIGN™ <i>Salmonella</i> -Ab	104912	strips	2	192	184
VETSIGN™ <i>Salmonella</i> -Ab	104913	plate	10	960	920

## ASFV

### African Swine Fever Virus Antibody test

An indirect ELISA for the detection of specific antibodies to African Swine Fever Virus in porcine serum samples (domestic and wild boars).

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® ASFV-Ab confirmation	104907	strips	2	96	88

## CSFV

### Classical Swine Fever Virus Antibody test

An indirect ELISA for the detection of specific antibodies to Classical Swine Fever Virus in porcine serum and plasma samples (domestic and wild boars).

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® CSFV-Ab	113135	strips	2	192	184
SVANOVIR® CSFV-Ab	113136	plate	10	960	920

## PRV

### **Pseudorabies/Aujeszky Disease Virus gE Antibody Test**

A blocking ELISA for the detection of specific antibodies to Pseudorabies/Aujeszky Disease Virus gE in porcine serum and plasma samples.

<b>Product name</b>	<b>Article No.</b>	<b>Format</b>	<b>Plates</b>	<b>Tests</b>	<b>Samples</b>
SVANOVIR® PRV gE-Ab	104902	strips	2	192	184
SVANOVIR® PRV gE-Ab	104904	plate	10	960	920

### **Pseudorabies/Aujeszky Disease Virus gB Antibody Test**

A blocking ELISA for the detection of specific antibodies to Pseudorabies/Aujeszky Disease Virus gB in porcine serum and plasma samples.

<b>Product name</b>	<b>Article No.</b>	<b>Format</b>	<b>Plates</b>	<b>Tests</b>	<b>Samples</b>
SVANOVIR® PRV gB-Ab	104905	strips	2	192	184
SVANOVIR® PRV gB-Ab	104906	plate	10	960	920

## BRUCELLA

### **Brucella C-ELISA Antibody Test**

A competitive ELISA for the detection of specific antibodies in porcine, bovine, ovine and caprine serum samples.

<b>Product name</b>	<b>Article No.</b>	<b>Format</b>	<b>Plates</b>	<b>Tests</b>	<b>Samples</b>
SVANOVIR® <i>Brucella</i> -Ab C-ELISA	104893	strips	2	192	176
SVANOVIR® <i>Brucella</i> -Ab C-ELISA	104894	plate	10	960	880

*Tests - number of tests. Samples - number of samples, wells for kit controls excluded.*

## EHV

### Equine Herpesvirus type 1 and 4 Discriminating test

An indirect ELISA for the detection of and/or discrimination of specific antibodies to Equine Herpesvirus type 1 and type 4 in equine serum and plasma samples.

Product name	Article No.	Format	Plates	Tests	Samples
SVANOVIR® EHV1/EHV4-Ab	104901	strips	2	64	52



# AVIAN

## APV

### **Avian Pneumovirus Antibody Test**

A blocking ELISA for the detection of specific antibodies to Avian Pneumovirus in serum and egg yolk from various avian species.

<b>Product name</b>	<b>Article No.</b>	<b>Format</b>	<b>Plates</b>	<b>Tests</b>	<b>Samples</b>
SVANOVIR® APV-Ab	104873	strips	2	192	184

## MG

### **Mycoplasma gallisepticum Antibody Test**

A blocking ELISA for the detection of specific antibodies to *Mycoplasma gallisepticum* in serum and egg yolk from various avian species.

<b>Product name</b>	<b>Article No.</b>	<b>Format</b>	<b>Plates</b>	<b>Tests</b>	<b>Samples</b>
SVANOVIR® MG-Ab	104875	strips	2	192	184

## MS

### **Mycoplasma synoviae Antibody Test**

A blocking ELISA for the detection of specific antibodies to *Mycoplasma synoviae* in serum and egg yolk from various avian species.

<b>Product name</b>	<b>Article No.</b>	<b>Format</b>	<b>Plates</b>	<b>Tests</b>	<b>Samples</b>
SVANOVIR® MS-Ab	104876	strips	2	192	184

Tests - number of tests. Samples - number of samples, wells for kit controls excluded.



### **We test our tests**

The quality of our processes, products and services are continuously optimized and improved in accordance with our customer's needs and the high standards of veterinary diagnostic laboratories. All products are developed, evaluated and manufactured according to ISO 9001:2008 to ensure consistency and reliability of our tests. With a combination of long experience and professionalism, Boehringer Ingelheim Svanova offers customers an ongoing, personal support. Our competence is part of our product.